PRO powered by BRAUN

TRIZONE

WATERJET 500 - 1000 / WATERJET

Type 3727  
Toothbrush Type 3756

Deutsch	6
English	10
Français	14
Español	19
Português	23
Italiano	27
Nederlands	31
Dansk / Norsk	35
Svenska	39
Suomi	43
Ελληνικά	47

Ⓝ Ⓜ 00 800 27 28 64 63

Ⓝ 08 44 - 88 40 10

Ⓝ 0 800 731 1792

Ⓝ 1 800 509 448

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Ⓝ 0 800 14 592

Ⓝ 800 21172

Ⓝ 901 11 61 84

Ⓝ 808 20 00 33

Ⓝ 800 124 600

Ⓝ 0800-4 45 53 88

Ⓝ 70 15 00 13

Ⓝ 22 63 00 93

Ⓝ 020 - 21 33 21

Ⓝ 020 377 877

Ⓝ 800 801 3457

**Internet:**

[www.oralb.com](http://www.oralb.com)  
[www.braun.com](http://www.braun.com)  
[www.service.braun.com](http://www.service.braun.com)

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# English

Welcome to Oral-B!

Please read the usage instructions carefully before first use and save this manual for future reference.

## IMPORTANT

- The appliance must be mounted in a way that it cannot fall into water or drop to the floor. The plug should always be accessible.
- Periodically check the cord for damage. If cord is damaged, take the charging unit to an Oral-B Service Centre. A damaged or non-functioning unit should no longer be used. Do not modify or repair the product. This may cause fire, electric shock or injury.
- The toothbrush is not intended for use by children under age three. Children aged from 3 years and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge can use the toothbrush, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Cleaning and user maintenance shall not be performed by children.
- The oral irrigator is not intended for use by children under age 8. Children aged from 8 years and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge can use the product, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Oral irrigators shall be used only with liquids defined in the manufacturer's instructions.

- Use this product only for its intended use as described in this manual. Do not use attachments or chargers which are not recommended by the manufacturer.

## WARNING

- If the product is dropped, the brush head/irrigator jet should be replaced before the next use even if no damage is visible.
- Do not place the appliance in water or liquid or store where it can fall or be pulled into a tub or sink. Do not reach for it when fallen into water. Unplug immediately.
- Do not disassemble the product except when disposing of the battery. When taking out the battery from the toothbrush for disposal, use caution not to short the positive (+) and negative (-) terminals.
- When unplugging, always hold the power plug instead of the cord. Do not touch the power plug with wet hands. This can cause an electric shock.
- If you are undergoing treatment for any oral care condition, consult your dental professional prior to use.
- This appliance is a personal care device and is not intended for use on multiple patients in a dental practice or institution.

## Description

### Irrigator

- 1 Jet compartment
- 2 Irrigator jet
- 3 Jet release button
- 4 On/off switch
- 5 Water pressure control dial
- 6 Handle switch (water stop and go)
- 7 Irrigator handle
- 8 Water container (600 ml)
- 9 Cord storage

### Toothbrush (depending on model)

- 10 Charger base
- 11 Low Charge Indicator light
- 12 Charge Indicator light
- 13 Toothbrush handle (rechargeable)
- 14 On/off button
- 15 Brush head
- 16 Brush head compartment

## Specifications

For voltage specifications please refer to the bottom of the appliance.

## A Using your irrigator

Before using the irrigator first time, fill the water container (8) and allow the water to run through

completely by holding the irrigator handle (7) over a washbasin.

- Attach one of the irrigator jets (2) to the handle until it locks into place (a). Fill the water container with lukewarm water. If you want to add mouthwash into the water, please make sure you fill the water container with lukewarm water first, then add the mouthwash (b).
- Turn on the irrigator by pressing the on/off switch (4) (c). The handle switch should be in the lower position = water stop.
- Select the water pressure with the pressure control dial (5) (1 = gentle, 5 = strong) (d). We recommend that you start with setting 1, especially if your gums tend to bleed or if children are using the irrigator.
- Bend over a washbasin and guide the jet to your teeth and gums. Slide the handle switch (6) up to start the water stream (e). Keep your mouth slightly open to allow the water to run into the washbasin.

The irrigator jet offers two different functions depending on the position of the jet tip:

#### Single jet:

Jet tip is in position (g) to clean teeth and the interdental spaces and remove resistant food particles.

#### Multi jet:

Jet tip is in position (h) to massage the gums and keep them healthy.

First, turn on the water stream by sliding the handle switch down (to avoid unintended water spraying). The jet can easily be switched from one function to the other by pressing the jet tip in or out. Begin with the single jet. Guide the water stream to the teeth; especially to the interdental spaces. Then massage the gums with the multi jet.

Maximal operation time: 15 minutes

Cooling off time: 2 hours

In general we recommend to replace the irrigator jet every 3 months.

#### **After use**

Always empty the water container entirely. If your dentist has recommended the use of mouthwash with your irrigator, allow some clean water to run through the irrigator before you empty the container entirely to prevent any deposit build-up. Make sure you keep the irrigator switched on until the water stops running through the irrigator jet, to prevent bacteria from collecting in the standing water.

Press the on/off switch to turn off the irrigator.

To take off the irrigator jet, press the jet release button (3) (f). Always wipe the handle and the

irrigator jet dry; place the handle back onto its base and store the irrigator jet in the jet compartment (1).

#### **Cleaning recommendations**

Regularly clean the housing with a damp cloth. The water container and the jet compartment are removable and dishwasher-safe.

## **B Using your toothbrush (depending on model)**

#### **Connecting and charging**

Your toothbrush has a waterproof handle, is electrically safe and designed for use in the bathroom.

- Plug the appliance into an electrical outlet and place the toothbrush handle (13) on the charge base (10) (i).
- The green charge indicator light (12) flashes while the handle is being charged. Once it is fully charged, the light turns off. A full charge can take up to 22 hours and enables up to seven days of regular brushing (twice a day, 2 minutes).  
Note: After deep discharge, the green light might not flash immediately; it can take up to 10–15 minutes.
- If the rechargeable battery is low the red low charge indicator light (11) is flashing for some seconds when turning your toothbrush on/off.
- For everyday use, the toothbrush handle can be stored on the Charger base to maintain it at full power. Overcharging is impossible.
- To maintain the maximum capacity of the rechargeable battery, unplug the charger base and fully discharge the handle with regular use at least every 6 months.

#### **Brushing techniques**

Wet brush head and apply any kind of toothpaste. To avoid splashing, guide the brush head (15) to your teeth before switching on the appliance (j). When brushing your teeth with one of the **Oral-B oscillating-rotating brush heads** guide the brush head slowly from tooth to tooth, spending a few seconds on each tooth surface (n). When using the **Oral-B «TriZone» brush head** place the toothbrush bristles against the teeth at a slight angle towards the gumline. Apply light pressure and start brushing in back and forth motions, just like you would do with a manual toothbrush. With any brush head start brushing the outsides, then the insides and finally the chewing surfaces. Brush all four quadrants of your mouth equally. Do not press too hard, simply let the brush do all the work. You may also consult your dentist or dental hygienist about the right technique for you.

In the first days of using any electric toothbrush, your gums may bleed slightly. In general, bleeding should stop after a few days. Should it persist after 2 weeks,

please consult your dentist or dental hygienist. If you have sensitive teeth and gums, Oral-B recommends using the «Sensitive» mode (depending on model) (optional in combination with an Oral-B «Sensitive» brush head).

### Brushing modes (depending on model)

- «Daily Clean» – Standard mode for daily mouth cleaning
- «Sensitive» – Gentle, yet thorough cleaning for sensitive areas

How to switch between modes:  
Your toothbrush automatically starts in «Daily Clean» mode. To switch to the other mode, press the on/off button (14). If you wish to turn off your toothbrush from the «Daily Clean» mode, push and hold down the on/off button until the motor stops.

### Professional Timer

A short stuttering sound at 30-second intervals reminds you to brush all four quadrants of your mouth equally (k). A long stuttering sound indicates the end of the professionally recommended 2-minute brushing time. The elapsed brushing time is memorised even when the handle is briefly switched off during brushing. When pausing longer than 30 seconds the timer resets.

### Pressure sensor

To promote optimal brushing your toothbrush has a pressure control feature installed. If too much pressure is applied, the oscillating movement of the brush head will continue but the pulsation will stop. You will also hear a different sound while brushing.

### Brush heads

Oral-B offers you a variety of different brush heads that fit your Oral-B toothbrush handle. Our **oscillating-rotating brush heads** can be used for tooth-by-tooth precision cleaning

#### Oral-B CrossAction brush head

Our most advanced brush head. Angled bristles for a precise clean. Lifts and powers away plaque.

#### Oral-B FlossAction brush head

features micropulse bristles that allow for an outstanding plaque removal in the interdental areas.

#### Oral-B 3D White brush head

features a special polishing cup to naturally whiten your teeth. Please note that children under 12 years old should not use the Oral-B «3D White» brush head.

#### Oral-B Sensitive brush head

features an extra soft bristle structure that is gentle on teeth and gums.

#### Oral-B PowerTip brush head

offers a deep clean in hard to reach areas, e.g. around bridges, crowns or implants.

Our **sweeping «Trizone» brush head** can be used with all common brushing techniques.

#### Oral-B TriZone brush head

features triple zone cleaning action for outstanding plaque removal, even in between teeth.

Oral-B brush heads feature light blue INDICATOR® bristles to help you monitor brush head replacement need. With thorough brushing, twice a day for two minutes, the blue color will fade halfway approximately within 3 months, indicating the need to replace your brush head. If the bristles splay before the color recedes, you may exert too much pressure on teeth and gums (m).

We do not recommend using the Oral-B FlossAction or Oral-B 3D White brush head with braces. You may use the Oral-B Ortho brush head, specifically designed to clean around braces and wires.

### Cleaning recommendations

After brushing, rinse your brush head carefully under running water for some seconds with the handle switched on. Switch off the handle and remove the brush head. Clean both parts separately under running water (l); then wipe them dry before reassembling the toothbrush. Put the toothbrush handle back on the charger base. The charger base should be cleaned with a damp cloth only.

Subject to change without notice.

### Environmental Notice

Product contains batteries and/or recyclable electric waste. For environment protection do not dispose in household waste, but for recycling take to electric waste collection points provided in your country.

## Guarantee

We grant 2 years guarantee on the product commencing on the date of purchase. Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose.

This guarantee extends to every country where this appliance is supplied by Braun/Oral-B or its appointed distributor.

This guarantee does not cover: damage due to improper use, normal wear or use, especially regarding the brush heads, as well as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorised persons and if original Braun parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Oral-B Braun Customer Service Centre.

### For UK only:

This guarantee in no way affects your rights under statutory law.

## Warranty replacement brush heads

Oral-B warranty will be voided if electric rechargeable handle damage is found to be attributed to the use of non-Oral-B replacement brush heads.

Oral-B does not recommend the use of non Oral-B replacement brush heads.

- Oral-B has no control over the quality of non-Oral-B replacement brush heads. Therefore, we cannot ensure the cleaning performance of non-Oral-B replacement brush heads, as communicated with the electric rechargeable handle at time of initial purchase.
- Oral-B cannot ensure a good fit of non-Oral-B replacement brush heads.
- Oral-B cannot predict the long-term effect of non-Oral-B replacement brush heads on handle wear.

All Oral-B replacement brush heads carry the Oral-B logo and meet Oral-B's high quality standards.

Oral-B does not sell replacement brush heads or handle parts under any other brand name.

## 30-Days Money Back Guarantee

Take the Oral-B 30 day Challenge: Satisfaction guaranteed or your money back!

Full details at [www.oralb.com](http://www.oralb.com) or call:

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